

Diamond Burs, Pack of 5


Flat End Cylindrical

- Useful for smoothing bone
- Cavity preparation
- Crown and bridge preparation

Item Number	Description	Bur
DTP7127C	Coarse	
DTP835-014	Coarse Feline	
DTP837-018	Coarse Canine	





Tight Occlusion

- Controlled contouring and finishing
- Smoothing bone

Item Number	Description	Bur
DTP7611TF	Trimming and Finishing/TF 12 Flutes	




Flame

- Cavity preparation
- Crown and bridge preparation

Item Number	Description	Bur
DTP5128M	Medium	
DTP5128C	Coarse	
DTP5128F	Fine	
DTP5128VF	Very Fine	




Round End Taper

- Crown and bridge preparation
- Round end designed to widen canal wall without the risk of gouging

Item Number	Description	Bur
DTP1610M	Medium	
DTP1610C	Coarse	
DTP1610F	Fine	


Egg

- Preparation of palatal areas
- Occlusal reduction
- Cavity preparation
- Crown and bridge preparation

Item Number	Description	Bur
DTP1908C	Coarse	
DTP1908F	Fine	
DTP390VF	Very Fine	

Football

- Preparation of palatal areas
- Occlusal reduction
- Cavity preparation
- Crown and bridge preparation

Item Number	Description	Bur
DTP1923M	Medium	

*All bur photos are enlarged to show detail.