



Friction Grip Burs, 16mm, Pack of 5





Straight Fissure Cross-Cut

All purpose burs that can be used to section teeth, gain endodontic access, and serve as a general-purpose bur.

Item Number	Cat. No.	Burs
DTP10021	#557	
DTP10019	#558	




Inverted Cone

Utilized for the undercutting cavity preparation and to produce a mechanical interlock.

Item Number	Cat. No.	Burs
DTP10033	#33 1/2	
DTP10034	#34	
DTP10035	#35	
DTP10037	#37	






Pear

Used for root canal access and for undercutting dentin.

Item Number	Cat. No.	Burs
DTP10017	#330	
DTP10018	#331	
DTP10020	#332	








Tapered Fissure Cross-Cut

All-purpose bur that can be used to section teeth, gain endodontic access, and serve as a general-purpose bur.

Item Number	Cat. No.	Burs
DTP10025	#699	
DTP10026	#700	
DTP10028	#701	
DTP10030	#701L	
DTP10032	#702	

Round

Used to gain endodontic access, to prepare cavities and to mark the cutting depth in crown preparation.

Item Number	Cat. No.	Burs
DTP10046	#1/2	
DTP10047	#1	
DTP10048	#2	
DTP10050	#3	
DTP10052	#4	
DTP10054	#5	
DTP10056	#6	

*All bur photos are enlarged to show detail.