




HP Burs, 40mm, Pack of 5






Tapered Fissure Cross-Cut

All-purpose bur that can be used to section teeth, gain endodontic access, and serve as a general-purpose bur.

Item Number	Cat. No.	Burs
DTP10070	#700	
DTP10071	#701	
DTP10072	#702	



Round

Used to gain endodontic access, to prepare cavities and to mark the cutting depth in crown preparation.

Item Number	Cat. No.	Burs
DTP10073	#1/2	
DTP10074	#2	
DTP10075	#4	
DTP10076	#6	
DTP10077	#8	

Inverted Cone

Utilized for the undercutting cavity preparation and to produce a mechanical interlock.

Item Number	Cat. No.	Burs
DTP10033	#33 1/2	
DTP10035	#35	

Recommended DentalAire Handpieces



Low-Speed Handpiece
DA500A



Low-speed Lube Free Handpiece
DA700A

*All bur photos are enlarged to show detail.