




## RA Latch Type Burs, Pack of 5






### Tapered Fissure Cross-Cut

All-purpose bur that can be used to section teeth, gain endodontic access, and serve as a general-purpose bur.

Item Number	Cat. No.	Burs
<b>DTP10036</b>	#700	
<b>DTP10038</b>	#701	
<b>DTP10042</b>	#702	

### Round

Used to gain endodontic access, to prepare cavities and to mark the cutting depth in crown preparation.

Item Number	Cat. No.	Burs
<b>DTP10059</b>	#1/2	
<b>DTP10060</b>	#2	
<b>DTP10064</b>	#4	
<b>DTP10066</b>	#6	
<b>DTP10068</b>	#8	





### Latch Type Burs


- Use with **low speed contra-angle latch type** handpieces
- Easily hooks to latch on low-speed handpiece

### Inverted Cone

Utilized for the undercutting cavity preparation and to produce a mechanical interlock.

Item Number	Cat. No.	Burs
<b>DTP10040</b>	#33 1/2	
<b>DTP10067</b>	#35	

### Recommended DentalAire Handpieces



**Contra Angle,  
Latch Type**  
DA553UL



**E-type Contra Angle,  
Latch Type**  
DA2500ECA

\*All bur photos are enlarged to show detail.