





Surgical Length Burs, 25mm, Pack of 5



Round

Used to gain endodontic access, to prepare cavities and to mark the cutting depth in crown preparation.

Item Number	Cat. No.	Burs
DTPFGSL2	FC. #2	
DTPFGSL4	FC. #4	
DTPFGSL6	FC. #6	
DTPFGSL8	FC. #8	



Straight Fissure Cross-Cut

All purpose burs that can be used to section teeth, gain endodontic access, and serve as a general-purpose bur.

Item Number	Cat. No.	Burs
DTPFGSL557	FC. #557	
DTPFGSL558	FC. #558	

Tapered Fissure Cross-Cut

All-purpose bur that can be used to section teeth, gain endodontic access, and serve as a general-purpose bur.

Item Number	Cat. No.	Burs
DTPFGSL701	FC. #701	
DTPFGSL702	FC. #702	

Recommended DentalAire Handpieces



High-Speed Swivel LED Handpiece
DA400LEDMH



High-speed LED Handpiece
DA400LED



High-Speed Fiber-optic Swivel Handpiece
DA400AFOPB-SW



High-speed Fiber-optic Handpiece
DA400FOPB

*All bur photos are enlarged to show detail.



(800)866-6881



info@dentalaire.com



www.dentalaire.com