

Prestige or Ultimate Dental Station?

SOLUTIONS FOR EVERY CLINIC.

The **Prestige and Ultimate** mobile dental stations help improve veterinary dental procedures for clinics and ensure that each patient receives the highest quality of care they deserve.



Unlimited Customer and Technical Support



Included DentaLaire Dental Vet Package

PRESTIGE DENTAL STATION



ULTIMATE DENTAL STATION



Design	Asepsis, Modern with Pneumatic Lift <ul style="list-style-type: none"> Helps promote a clean working station and eliminate harmful bacteria 	Traditional with Manual Height Adjustment <ul style="list-style-type: none"> Ergonomic
Handpieces	LED Piezo Scaler <ul style="list-style-type: none"> LED Piezo Scaler, LED Swivel High-speed, Low-speed, Air/Water Syringe 	Piezo Scaler, with No Light <ul style="list-style-type: none"> Piezo Scaler (with no light), LED Swivel High-speed, Low-speed, Air/Water Syringe
Water Bottle System	2 - 2 Liter Water Bottles <ul style="list-style-type: none"> Reliable water source during dental procedures and cleanings 	2 - 2 Liter Water Bottles <ul style="list-style-type: none"> Reliable water source during dental procedures and cleanings
Operation Controls	Side Located Recessed <ul style="list-style-type: none"> Enhanced ergonomics, improved accessibility, and maximized workspace 	Located on Front Panel <ul style="list-style-type: none"> Easy to access controls and user-friendly
Storage	Internal <ul style="list-style-type: none"> Large storage tray located under lid 	No Storage <ul style="list-style-type: none"> Large workspace area that includes removable stainless steel tray
Service	Enhanced Serviceability <ul style="list-style-type: none"> All service parts are easily accessible under storage tray in lid 	Service Parts Under Top Countertop <ul style="list-style-type: none"> Access service parts under top countertop
Compressor	Pick your option <ul style="list-style-type: none"> Oil, oil-less, and no compressor option 	Pick your option <ul style="list-style-type: none"> Oil, oil-less, and no compressor option
Warranty	2 year warranty* <ul style="list-style-type: none"> Excludes handpieces and accessories 	1 year warranty* <ul style="list-style-type: none"> Excludes handpieces and accessories

